


Supplementary Table 1. SYBR green assay sequences

Gene	Forward primer	Reverse primer
<i>PPARGC1A</i>	AGCCTCTTTGCCCAGATCTT	GGCAATCCGTCTTCATCCAC
<i>36B4</i>	GTGATGTGCAGCTGATCAAGACT	GATGACCAGCCCAAAGGAGA
<i>BAX</i>	AGCTGCAGAGGATGATTGCC	GCGTCCCAAAGTAGGAGAGG
<i>B2M</i>	GTGATGTGCAGCTGATCAAGACT	GATGACCAGCCCAAAGGAGA

Supplementary Table 2. Taqman assay ID's

Gene	Assay ID
<i>MT-RNR2</i>	HS02596860s1
<i>Humanin</i>	Custom designed Taqman small RNA assay against humanin sequence https://www.thermofisher.com/order/catalog/product/4398987?SID=srch-srp-4398987#/4398987?SID=srch-srp-4398987
<i>RNU44</i>	001094
<i>RNU48</i>	001006



Supplementary figure 1. The custom TaqMan® small RNA assays designed against the humanin sequence is not amplified in cells that do not have any mitochondrial DNA (143B Rho 0 cells) **(a)**, and humanin mRNA levels in large (>200 nucleotides) and small (<200 nucleotides) mRNA muscle fractions at rest before initiation of exercise training **(b)**. Absolute $\Delta\Delta$ CT data for *Humanin* mRNA in small **(c)** and large fraction **(d)**. *** $p < 0.001$ vs small fraction for student t-test.